

(19) World Intellectual Property Organization  
International Bureau



(43) International Publication Date  
21 March 2002 (21.03.2002)

PCT

(10) International Publication Number  
**WO 02/22828 A1**

(51) International Patent Classification<sup>7</sup>: - - G12N 15/54,  
15/74, 1/21, 9/12, C12P 13/04, 13/08, C12Q 1/68

(21) International Application Number: PCT/EP01/10211

(22) International Filing Date:  
5 September 2001 (05.09.2001)

(25) Filing Language: English

(26) Publication Language: English

(30) Priority Data:  
100 44 912.3 12 September 2000 (12.09.2000) DE  
101 20 095.1 25 April 2001 (25.04.2001) DE

(71) Applicant: DEGUSSA AG [DE/DE]; Bennigsenplatz 1,  
40474 Düsseldorf (DE).

(72) Inventors: BATHE, Brigitte; Twieten 1, 33154 Salzkot-  
ten (DE). HANS, Stephan; Wilhelmstrasse 7, 49078 Os-  
nabrück (DE). FARWICK, Mike; Gustav-Adolf-Strasse  
11, 33615 Bielefeld (DE). HERMANN, Thomas; Zirkon-  
strasse 8, 33739 Bielefeld (DE).

(81) Designated States (national): AE, AG, AL, AM, AT, AU,  
AZ, BA, BB, BG, BR, BY, BZ, CA, CH, CN, CO, CR, CU,

CZ, DE, DK, DM, DZ, EC, EE, ES, FI, GB, GD, GE, GH,  
GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC,  
LK, LR, LS, LT, LU, LV, MA, MD, MG, MK, MN, MW,  
MX, MZ, NO, NZ, PH, PL, PT, RO, RU, SD, SE, SG, SI,  
SK, SL, TJ, TM, TR, TT, TZ, UA, UG, UZ, VN, YU, ZA,  
ZW.

(84) Designated States (regional): ARIPO patent (GH, GM,  
KE, LS, MW, MZ, SD, SL, SZ, TZ, UG, ZW), Eurasian  
patent (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM), European  
patent (AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE,  
IT, LU, MC, NL, PT, SE, TR), OAPI patent (BF, BJ, CF,  
CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD,  
TG).

**Published:**

- with international search report
- before the expiration of the time limit for amending the  
claims and to be republished in the event of receipt of  
amendments
- with (an) indication(s) in relation to deposited biological  
material furnished under Rule 13bis separately from the  
description

For two-letter codes and other abbreviations, refer to the "Guid-  
ance Notes on Codes and Abbreviations" appearing at the begin-  
ning of each regular issue of the PCT Gazette.

WO 02/22828 A1

(54) Title: ISOLATION AND SEQUENCING OF THE PKNB GENE OF C. GLUTAMICUM

(57) Abstract: The invention relates to the protein kinase B of C. glutamicum and a fermentation process for the preparation of L-amino acids using corynebacteria in which at least the pknB gene is amplified, and to the use, as hybridization probes, of polynucleotides containing the sequences according to the invention.

## ISOLATION AND SEQUENCING OF THE PKNB GENE OF C. GLUTAMICUM

## Field of the Invention

The invention provides nucleotide sequences from corynebacteria coding for the pknB gene and a fermentation process for the preparation of amino acids using bacteria in which the endogeny pknB gene is amplified.

## State of the Art

L-Amino acids, especially L-lysine, are used in human medicine, in the pharmaceutical industry, in the food industry and very especially in animal nutrition.

It is known that amino acids are prepared by the fermentation of strains of corynebacteria, especially *Corynebacterium glutamicum*. Because of their great importance, attempts are constantly being made to improve the preparative processes. Improvements to the processes may relate to measures involving the fermentation technology, e.g. stirring and oxygen supply, or the composition of the nutrient media, e.g. the sugar concentration during fermentation, or the work-up to the product form, e.g. by ion exchange chromatography, or the intrinsic productivity characteristics of the microorganism itself.

The productivity characteristics of these microorganisms are improved by using methods of mutagenesis, selection and mutant choice to give strains which are resistant to antimetabolites or auxotrophic for metabolites important in regulation, and produce amino acids.

Methods of recombinant DNA technology have also been used for some years to improve L-amino acid-producing strains of *Corynebacterium* by amplifying individual amino acid biosynthesis genes and studying the effect on amino acid production.

## Object of the Invention

The object which the inventors set themselves was to provide novel measures for improving the preparation of amino acids by fermentation.

## Summary of the Invention

When L-amino acids or amino acids are mentioned hereafter, they are understood as meaning one or more amino acids, including their salts, selected from the group comprising L-asparagine, L-threonine, L-serine, L-glutamate, L-glycine, L-alanine, L-cysteine, L-valine, L-methionine, L-isoleucine, L-leucine, L-tyrosine, L-phenylalanine, histidine, L-lysine, L-tryptophan and L-arginine. L-lysine is particularly preferred.

When L-lysine or lysine is mentioned hereafter, it is understood as meaning not only the bases but also the salts, e.g. lysine monohydrochloride or lysine sulfate.

The invention provides an isolated polynucleotide from corynebacteria which contains a polynucleotide sequence coding for the pknB gene and is selected from the group comprising:

- a) a polynucleotide which is at least 70% identical to a polynucleotide coding for a polypeptide containing the amino acid sequence of SEQ ID No. 2,
- b) a polynucleotide coding for a polypeptide containing an amino acid sequence which is at least 70% identical to the amino acid sequence of SEQ ID No. 2,
- c) a polynucleotide which is complementary to the polynucleotides of a) or b), and
- d) a polynucleotide containing at least 15 consecutive nucleotides of the polynucleotide sequence of a), b) or c),

the polypeptide preferably having the activity of protein kinase B.

The invention also provides the above-mentioned polynucleotide, which is preferably a replicatable DNA containing:

- (i) the nucleotide sequence shown in SEQ ID No. 1, or
- (ii) at least one sequence corresponding to sequence (i) within the degeneracy of the genetic code, or
- (iii) at least one sequence which hybridizes with the sequence complementary to sequence (i) or (ii), and optionally
- (iv) neutral sense mutations in (i).

The invention also provides:

a replicatable polynucleotide, especially DNA, containing the nucleotide sequence as shown in SEQ ID No. 1,

a polynucleotide coding for a polypeptide containing the amino acid sequence as shown in SEQ ID No. 2,

a vector containing the polynucleotide according to the invention, especially a shuttle vector or plasmid vector, and

corynebacteria which contain the vector or in which the endogeny *pknB* gene is amplified.

The invention also provides polynucleotides consisting substantially of a polynucleotide sequence which are obtainable by screening, by means of hybridization, of an appropriate gene library of a corynebacterium, containing the complete gene or parts thereof, with a probe containing the sequence of the polynucleotide of the invention

according to SEQ ID No. 1 or a fragment thereof, and by isolation of said polynucleotide sequence.

#### Detailed Description of the Invention

As hybridization probes for RNA, cDNA and DNA, polynucleotides containing the sequences according to the invention are suitable for isolating the full length of nucleic acids, or polynucleotides or genes, coding for protein kinase B, or for isolating nucleic acids, or polynucleotides or genes, whose sequence exhibits a high degree of similarity to the sequence of the pknB gene. They are also suitable for incorporation into so-called arrays, micro-arrays or DNA chips for detecting and determining the corresponding polynucleotides.

Polynucleotides containing the sequences according to the invention are further suitable as primers for the preparation, by the polymerase chain reaction (PCR), of DNA of genes coding for protein kinase B.

Such oligonucleotides serving as probes or primers contain at least 25, 26, 27, 28, 29 or 30, preferably at least 20, 21, 22, 23 or 24 and very particularly preferably at least 15, 16, 17, 18 or 19 consecutive nucleotides. Oligonucleotides with a length of at least 31, 32, 33, 34, 35, 36, 37, 38, 39 or 40 or at least 41, 42, 43, 44, 45, 46, 47, 48, 49 or 50 nucleotides are also suitable. Oligonucleotides with a length of at least 100, 150, 200, 250 or 300 nucleotides may also be suitable.

"Isolated" means separated from its natural environment.

"Polynucleotide" refers in general to polyribonucleotides and polydeoxyribonucleotides, it being possible for the RNAs or DNAs in question to be unmodified or modified.

The polynucleotides according to the invention include a polynucleotide according to SEQ ID No. 1 or a fragment

prepared therefrom, as well as polynucleotides which are in particular at least 70% to 80%, preferably at least 81% to 85%, particularly preferably at least 86% to 90% and very particularly preferably at least 91%, 93%, 95%, 97% or 99% identical to the polynucleotide according to SEQ-ID No. 1 or a fragment prepared therefrom.

"Polypeptides" are understood as meaning peptides or proteins containing two or more amino acids bonded via peptide links.

The polypeptides according to the invention include a polypeptide according to SEQ ID No. 2, especially those with the biological activity of protein kinase B and also those which are at least 70% to 80%, preferably at least 81% to 85%, particularly preferably at least 86% to 90% and very particularly preferably at least 91%, 93%, 95%, 97% or 99% identical to the polypeptide according to SEQ ID No. 2, and have said activity.

The invention further relates to a fermentation process for the preparation of amino acids selected from the group comprising L-asparagine, L-threonine, L-serine, L-glutamate, L-glycine, L-alanine, L-cysteine, L-valine, L-methionine, L-isoleucine, L-leucine, L-tyrosine, L-phenylalanine, L-histidine, L-lysine, L-tryptophan and L-arginine, using corynebacteria which, in particular, already produce amino acids and in which the nucleotide sequences coding for the *pknB* gene are amplified and, in particular, overexpressed.

In this context the term "enhancement" describes the increase in the intracellular activity, in a microorganism, of one or more enzymes which are coded for by the appropriate DNA, for example by increasing the copy number of the gene(s) or allele(s), using a strong promoter or using a gene or allele coding for an appropriate enzyme

with a high activity, and optionally combining these measures..

By enhancement measures, in particular over-expression, the activity or concentration of the corresponding protein is in general increased by at least 10%, 25%, 50%, 75%, 100%, 150%, 200%, 300%, 400% or 500%, up to a maximum of 1000% or 2000%, based on the starting microorganism.

The microorganisms provided by the present invention can produce L-amino acids from glucose, sucrose, lactose, fructose, maltose, molasses, starch or cellulose or from glycerol and ethanol. Said microorganisms can be representatives of corynebacteria, especially of the genus *Corynebacterium*. The species *Corynebacterium glutamicum* may be mentioned in particular in the genus *Corynebacterium*, being known to those skilled in the art for its ability to produce L-amino acids.

The following known wild-type strains:

- Corynebacterium glutamicum* ATCC13032
- Corynebacterium acetoglutamicum* ATCC15806
- Corynebacterium acetoacidophilum* ATCC13870
- Corynebacterium thermoaminogenes* FERM BP-1539
- Corynebacterium melassecola* ATCC17965
- Brevibacterium flavum* ATCC14067
- Brevibacterium lactofermentum* ATCC13869 and
- Brevibacterium divaricatum* ATCC14020

and L-amino acid-producing mutants or strains prepared therefrom, are particularly suitable strains of the genus *Corynebacterium*, especially of the species *Corynebacterium glutamicum* (*C. glutamicum*).

The novel *pknB* gene of *C. glutamicum* coding for the enzyme protein kinase B (EC 2.7.1.37) has been isolated.

The first step in isolating the *pknB* gene or other genes of *C. glutamicum* is to construct a gene library of this microorganism in *Escherichia coli* (*E. coli*). The construction of gene libraries is documented in generally well-known textbooks and manuals. Examples which may be mentioned are the textbook by Winnacker entitled *From Genes to Clones, Introduction to Gene Technology* (Verlag Chemie, Weinheim, Germany, 1990) or the manual by Sambrook et al. entitled *Molecular Cloning, A Laboratory Manual* (Cold Spring Harbor Laboratory Press, 1989). A very well-known gene library is that of the *E. coli* K-12 strain W3110, which was constructed by Kohara et al. (*Cell* 50, 495-508 (1987)) in  $\lambda$  vectors. Bathe et al. (*Molecular and General Genetics* 252, 255-265, 1996) describe a gene library of *C. glutamicum* ATCC13032, which was constructed using cosmid vector SuperCos I (Wahl et al., 1987, *Proceedings of the National Academy of Sciences USA* 84, 2160-2164) in the *E. coli* K-12 strain NM554 (Raleigh et al., 1988, *Nucleic Acids Research* 16, 1563-1575).

Börmann et al. (*Molecular Microbiology* 6(3), 317-326 (1992)) in turn describe a gene library of *C. glutamicum* ATCC13032 using cosmid pH79 (Hohn and Collins, *Gene* 11, 291-298 (1980)).

A gene library of *C. glutamicum* in *E. coli* can also be constructed using plasmids like pBR322 (Bolivar, *Life Sciences* 25, 807-818 (1979)) or pUC9 (Vieira et al., 1982, *Gene* 19, 259-268). Restriction- and recombination-defective *E. coli* strains are particularly suitable as hosts, an example being the strain DH5 $\alpha$ mcr, which has been described by Grant et al. (*Proceedings of the National Academy of Sciences USA* 87 (1990) 4645-4649). The long DNA fragments cloned with the aid of cosmids can then in turn be subcloned into common vectors suitable for sequencing, and subsequently sequenced, e.g. as described by Sanger et

al. (Proceedings of the National Academy of Sciences of the United States of America 74, 5463-5467, 1977).

The DNA sequences obtained can then be examined with known algorithms or sequence analysis programs, e.g. that of Staden (Nucleic Acids Research 14, 217-232 (1986)), that of Marck (Nucleic Acids Research 16, 1829-1836 (1988)) or the GCG program of Butler (Methods of Biochemical Analysis 39, 74-97 (1998)).

The novel DNA sequence of *C. glutamicum* coding for the *pknB* gene was found and, as SEQ ID No. 1, forms part of the present invention. Furthermore, the amino acid sequence of the corresponding protein was derived from said DNA sequence by the methods described above. The resulting amino acid sequence of the *pknB* gene product is shown in SEQ ID No. 2.

Coding DNA sequences which result from SEQ ID No. 1 due to the degeneracy of the genetic code also form part of the invention. Likewise, DNA sequences which hybridize with SEQ ID No. 1 or portions of SEQ ID No. 1 form part of the invention. Furthermore, conservative amino acid exchanges, e.g. the exchange of glycine for alanine or of aspartic acid for glutamic acid in proteins, are known to those skilled in the art as "sense mutations", which do not cause a fundamental change in the activity of the protein, i.e. they are neutral. It is also known that changes at the N- and/or C-terminus of a protein do not substantially impair its function or may even stabilize it. Those skilled in the art will find information on this subject in Ben-Bassat et al. (Journal of Bacteriology 169, 751-757 (1987)), O'Regan et al. (Gene 77, 237-251 (1989)), Sahin-Toth et al. (Protein Sciences 3, 240-247 (1994)) and Hochuli et al. (Bio/Technology 6, 1321-1325 (1988)), inter alia, and in well-known textbooks on genetics and molecular biology. Amino acid sequences which correspondingly result from SEQ ID No. 2 also form part of the invention.

Likewise, DNA sequences which hybridize with SEQ ID No. 1 or portions of SEQ ID No. 1 form part of the invention. Finally, DNA sequences which are prepared by the polymerase chain reaction (PCR) using primers resulting from SEQ ID No. 1 form part of the invention. Such oligonucleotides typically have a length of at least 15 nucleotides.

Those skilled in the art will find instructions on the identification of DNA sequences by means of hybridization in the manual entitled "The DIG System User's Guide for Filter Hybridization" from Boehringer Mannheim GmbH (Mannheim, Germany, 1993) and in Liebl et al. (International Journal of Systematic Bacteriology (1991) 41, 255-260), inter alia. Hybridization takes place under stringent conditions; in other words, only hybrids for which the probe and the target sequence, i.e. the polynucleotides treated with the probe, are at least 70% identical are formed. It is known that the stringency of hybridization, including the washing steps, is influenced or determined by varying the buffer composition, the temperature and the salt concentration. The hybridization reaction is preferably carried out under relatively low stringency compared with the washing steps (Hybaid Hybridisation Guide, Hybaid Limited, Teddington, UK, 1996).

The hybridization reaction can be carried out for example using a 5x SSC buffer at a temperature of approx. 50°C - 68°C, it also being possible for probes to hybridize with polynucleotides which are less than 70% identical to the sequence of the probe. Such hybrids are less stable and are removed by washing under stringent conditions. This can be achieved for example by lowering the salt concentration to 2x SSC and subsequently to 0.5x SSC if necessary (The DIG System User's Guide for Filter Hybridization, Boehringer Mannheim, Mannheim, Germany, 1995), the temperature being adjusted to approx. 50°C - 68°C. It is possible to lower the salt concentration to

0.1x SSC if necessary. By raising the hybridization temperature in approx. 1 - 2°C steps from 50°C to 68°C, it is possible to isolate polynucleotide fragments which are e.g. at least 70%, at least 80% or at least 90% to 95% identical to the sequence of the probe used. Further instructions on hybridization are commercially available in the form of kits (e.g. DIG Easy Hyb from Roche Diagnostics GmbH, Mannheim, Germany, Catalog No. 1603558).

Those skilled in the art will find instructions on the amplification of DNA sequences with the aid of the polymerase chain reaction (PCR) in the manual by Gait entitled Oligonucleotide Synthesis: A Practical Approach (IRL Press, Oxford, UK, 1984) and in Newton and Graham: PCR (Spektrum Akademischer Verlag, Heidelberg, Germany, 1994), inter alia.

It has been found that, after overexpression of the pknB gene, the production of amino acids by corynebacteria is improved.

Overexpression can be achieved by increasing the copy number of the appropriate genes or mutating the promoter and regulatory region or the ribosome binding site located upstream from the structural gene. Expression cassettes incorporated upstream from the structural gene work in the same way. Inducible promoters additionally make it possible to increase the expression in the course of the production of amino acid by fermentation. Measures for prolonging the life of the mRNA also improve the expression. Furthermore, the enzyme activity is also enhanced by preventing the degradation of the enzyme protein. The genes or gene constructs can either be located in plasmids of variable copy number or integrated and amplified in the chromosome. Alternatively, it is also possible to achieve overexpression of the genes in question by changing the composition of the media and the culture technique.

Those skilled in the art will find relevant instructions in Martin et al. (Bio/Technology 5, 137-146 (1987)), Guerrero et al. (Gene 138, 35-41 (1994)), Tsuchiya and Morinaga (Bio/Technology 6, 428-430 (1988)), Eikmanns et al. (Gene 102, 93-98 (1991)), EP 0 472 869, US 4,601,893, Schwarzer and Pühler (Bio/Technology 9, 84-87 (1991)), Reinscheid et al. (Applied and Environmental Microbiology 60, 126-132 (1994)), LaBarre et al. (Journal of Bacteriology 175, 1001-1007 (1993)), WO 96/15246, Malumbres et al. (Gene 134, 15-24 (1993)), JP-A-10-229891, Jensen and Hammer (Biotechnology and Bioengineering 58, 191-195 (1998)) and Makrides (Microbiological Reviews 60, 512-538 (1996)), inter alia, and in well-known textbooks on genetics and molecular biology.

For amplification, the *pknB* gene according to the invention has been overexpressed for example with the aid of episomal plasmids. Suitable plasmids are those which are replicated in corynebacteria. Numerous known plasmid vectors, e.g. pZ1 (Menkel et al., Applied and Environmental Microbiology (1989) 64, 549-554), pEKEx1 (Eikmanns et al., Gene 102, 93-98 (1991)) or pHS2-1 (Sonnen et al., Gene 107, 69-74 (1991)), are based on cryptic plasmids pHM1519, pBL1 or pGA1. Other plasmid vectors, e.g. those based on pCG4 (US-A-4,489,160), pNG2 (Serwold-Davis et al., FEMS Microbiology Letters 66, 119-124 (1990)) or pAG1 (US-A-5,158,891), can be used in the same way.

Other suitable plasmid vectors are those which make it possible to use the gene amplification process by integration into the chromosome, as described for example by Reinscheid et al. (Applied and Environmental Microbiology 60, 126-132 (1994)) for the duplication or amplification of the *hom-thrB* operon. In this method the complete gene is cloned into a plasmid vector which can replicate in a host (typically *E. coli*), but not in *C. glutamicum*. Examples of suitable vectors are pSUP301

(Simon et al., Bio/Technology 1, 784-791 (1983)), pK18mob or pK19mob (Schäfer et al., Gene 145, 69-73 (1994)), pGEM-T (Promega Corporation, Madison, WI, USA), pCR2.1-TOPO (Shuman (1994), Journal of Biological Chemistry 269, 32678-84; US-A-5,487,993), pCR®Blunt (Invitrogen, Groningen, The Netherlands; Bernard et al., Journal of Molecular Biology 234, 534-541 (1993)), pEM1 (Schrumpf et al., 1991, Journal of Bacteriology 173, 4510-4516) or pBGS8 (Spratt et al., 1986, Gene 41, 337-342). The plasmid vector containing the gene to be amplified is then transferred to the desired strain of *C. glutamicum* by conjugation or transformation. The method of conjugation is described for example in Schäfer et al. (Applied and Environmental Microbiology 60, 756-759 (1994)). Methods of transformation are described for example in Thierbach et al. (Applied Microbiology and Biotechnology 29, 356-362 (1988)), Dunican and Shivnan (Bio/Technology 7, 1067-1070 (1989)) and Tauch et al. (FEMS Microbiological Letters 123, 343-347 (1994)). After homologous recombination by means of a crossover event, the resulting strain contains at least two copies of the gene in question.

It has also been found that amino acid exchanges in the section between position 581 and position 587 of the amino acid sequence of protein kinase B, shown in SEQ ID No. 2, improve the production of amino acids, especially lysine, by corynebacteria.

Preferably, L-proline in position 584 is exchanged for any other proteogenic amino acid except L-proline, preferably for L-serine or L-threonine and very particularly preferably for L-serine.

SEQ ID No. 3 shows the base sequence of the pknB-1547 allele contained in the strain DM1547. The pknB-1547 allele codes for a protein whose amino acid sequence is shown in SEQ ID No. 4. The protein contains L-serine in position 584. The DNA sequence of the pknB-1547 allele

(SEQ ID No. 3) contains the base thymine in place of the base cytosine contained in the pknB wild-type gene (SEQ ID No. 1) in position 2343.

----- Mutagenesis can be carried out by conventional methods using mutagenic substances such as N-methyl-N'-nitro-N-nitrosoguanidine or ultraviolet light. Mutagenesis can also be carried out using in vitro methods such as treatment with hydroxylamine (Miller, J.H.: A Short Course in Bacterial Genetics. A Laboratory Manual and Handbook for Escherichia coli and Related Bacteria, Cold Spring Harbor Laboratory Press, Cold Spring Harbor, 1992). or mutagenic oligonucleotides (T.A. Brown: Gentechnologie für Einsteiger (Gene Technology for Beginners), Spektrum Akademischer Verlag, Heidelberg, 1993), or the polymerase chain reaction (PCR) as described in the manual by Newton and Graham (PCR, Spektrum Akademischer Verlag, Heidelberg, 1994).

The corresponding alleles or mutations are sequenced and introduced into the chromosome by the method of gene replacement, for example as described in Peters-Wendisch et al. (Microbiology 144, 915-927 (1998)) for the pyc gene of C. glutamicum, in Schäfer et al. (Gene 145, 69-73 (1994)) for the hom-thrB gene region of C. glutamicum or in Schäfer et al. (Journal of Bacteriology 176, 7309-7319 (1994)) for the cgl gene region of C. glutamicum. The corresponding alleles or the associated proteins can optionally be amplified in turn.

In addition it can be advantageous for the production of L-amino acids to amplify and, in particular, overexpress not only the pknB gene but also one or more enzymes of the particular biosynthetic pathway, the glycolysis, the anaplerosis, the citric acid cycle, the pentose phosphate cycle or the amino acid export, and optionally regulatory proteins.

Thus, for the production of L-amino acids, one or more genes selected from the following group can be amplified and, in particular, overexpressed in addition to amplification of the endogene *pknB* gene:

- the *dapA* gene coding for dihydrodipicolinate synthase (EP-B-0 197 335),
- the *gap* gene coding for glyceraldehyde 3-phosphate dehydrogenase (Eikmanns (1992), Journal of Bacteriology 174, 6076-6086),
- the *tpi* gene coding for triose phosphate isomerase (Eikmanns (1992), Journal of Bacteriology 174, 6076-6086),
- the *pgk* gene coding for 3-phosphoglycerate kinase (Eikmanns (1992), Journal of Bacteriology 174, 6076-6086),
- the *zwf* gene coding for glucose-6-phosphate dehydrogenase (JP-A-09224661),
- the *pyc* gene coding for pyruvate carboxylase (DE-A-198 31 609),
- the *lysC* gene coding for a feedback-resistant aspartate kinase (Accession no. P26512; EP-B-0387527; EP-A-0699759),
- the *lysE* gene coding for lysine export (DE-A-195 48 222),
- the *hom* gene coding for homoserine dehydrogenase (EP-A-0131171),
- the *ilvA* gene coding for threonine dehydratase (Möckel et al., Journal of Bacteriology (1992) 8065-8072) or the *ilvA(Fbr)* allele coding for a feedback-resistant

threonine dehydratase (Möckel et al. (1994), Molecular Microbiology 13, 833-842),

- the ilvBN gene coding for acetohydroxy acid synthase (EP-B-0356739),
- the ilvD gene coding for dihydroxy acid dehydratase (Sahm and Eggeling (1999), Applied and Environmental Microbiology 65, 1973-1979), or
- the zwal gene coding for the Zwal protein (DE 199 59 328.0, DSM13115).

In addition to amplification of the pknB gene, it can also be advantageous for the production of L-amino acids to attenuate one or more genes selected from the following group:

- the pck gene coding for phosphoenol pyruvate carboxykinase (DE 199 50 409.1, DSM13047),
- the pgi gene coding for glucose-6-phosphate isomerase (US 09/396,478, DSM12969),
- the poxB gene coding for pyruvate oxidase (DE 199 51 975.7, DSM13114), or
- the zwa2 gene coding for the Zwa2 protein (DE 199 59 327.2, DSM13113),

and, in particular, to reduce the expression.

In this context the term "attenuation" describes the reduction or switching-off of the intracellular activity, in a microorganism, of one or more enzymes (proteins) which are coded for by the appropriate DNA, for example by using a weak promoter or using a gene or allele coding for an appropriate enzyme with a low activity, or inactivating the appropriate gene or enzyme (protein), and optionally combining these measures.

By attenuation measures, the activity or concentration of the corresponding protein is in general reduced to 0 to 50%, 0 to 25%, 0 to 10% or 0 to 5% of the activity or concentration of the wild-type protein.

It can also be advantageous for the production of amino acids not only to overexpress the *pknB* gene but also to switch off unwanted secondary reactions (Nakayama: "Breeding of Amino Acid Producing Micro-organisms", in: Overproduction of Microbial Products, Krumphanzl, Sikyta, Vanek (eds.), Academic Press, London, UK, 1982).

The microorganisms prepared according to the invention are also provided by the invention and can be cultivated for the production of amino acids continuously or discontinuously by the batch process, the fed batch process or the repeated fed batch process. A summary of known cultivation methods is described in the textbook by Chmiel (Bioprozesstechnik 1. Einführung in die Bioverfahrenstechnik (Bioprocess Technology 1. Introduction to Bioengineering) (Gustav Fischer Verlag, Stuttgart, 1991)) or in the textbook by Storhas (Bioreaktoren und periphere Einrichtungen (Bioreactors and Peripheral Equipment) (Vieweg Verlag, Brunswick/Wiesbaden, 1994)).

The culture medium to be used must appropriately meet the demands of the particular strains. Descriptions of culture media for various microorganisms can be found in "Manual of Methods for General Bacteriology" of the American Society for Bacteriology (Washington DC, USA, 1981).

Carbon sources which can be used are sugars and carbohydrates, e.g. glucose, sucrose, lactose, fructose, maltose, molasses, starch and cellulose, oils and fats, e.g. soybean oil, sunflower oil, groundnut oil and coconut fat, fatty acids, e.g. palmitic acid, stearic acid and linoleic acid, alcohols, e.g. glycerol and thanol, and

organic acids, e.g. acetic acid. These substances can be used individually or as a mixture.

Nitrogen sources which can be used are organic nitrogen-containing compounds such as peptones, yeast extract, meat extract, malt extract, corn steep liquor, soybean flour and urea, or inorganic compounds such as ammonium sulfate, ammonium chloride, ammonium phosphate, ammonium carbonate and ammonium nitrate. The nitrogen sources can be used individually or as a mixture.

Phosphorus sources which can be used are phosphoric acid, potassium dihydrogenphosphate or dipotassium hydrogenphosphate or the corresponding sodium salts. The culture medium must also contain metal salts, e.g. magnesium sulfate or iron sulfate, which are necessary for growth. Finally, essential growth-promoting substances such as amino acids and vitamins can be used in addition to the substances mentioned above. Suitable precursors can also be added to the culture medium. Said feed materials can be added to the culture all at once or fed in appropriately during cultivation.

The pH of the culture is controlled by the appropriate use of basic compounds such as sodium hydroxide, potassium hydroxide, ammonia or aqueous ammonia, or acidic compounds such as phosphoric acid or sulfuric acid. Foaming can be controlled using antifoams such as fatty acid polyglycol esters. The stability of plasmids can be maintained by adding suitable selectively acting substances, e.g. antibiotics, to the medium. Aerobic conditions are maintained by introducing oxygen or oxygen-containing gaseous mixtures, e.g. air, into the culture. The temperature of the culture is normally 20°C to 45°C and preferably 25°C to 40°C. The culture is continued until the formation of the desired product has reached a maximum. This objective is normally achieved within 10 hours to 160 hours.

Methods of determining L-amino acids are known from the state of the art. They can be analyzed for example by ion exchange chromatography with subsequent ninhydrin derivatization, as described by Spackman et al. (Analytical Chemistry 30 (1958) 1190), or by reversed phase HPLC, as described by Lindroth et al. (Analytical Chemistry (1979) 51, 1167-1174).

A pure culture of the *Corynebacterium glutamicum* strain DM1547 was deposited as DSM13994 in the Deutsche Sammlung für Mikroorganismen und Zellkulturen (German Collection of Microorganisms and Cell Cultures (DSMZ), Brunswick, Germany) on 16 January 2001 under the terms of the Budapest Treaty.

The fermentation process according to the invention is used for the preparation of amino acids.

The present invention is illustrated in greater detail below by means of Examples.

The isolation of plasmid DNA from *Escherichia coli* and all the techniques of restriction, Klenow treatment and alkaline phosphatase treatment were carried out according to Sambrook et al. (Molecular Cloning. A Laboratory Manual (1989), Cold Spring Harbor Laboratory Press, Cold Spring Harbor, NY, USA). Methods of transforming *Escherichia coli* are also described in this manual.

The composition of common nutrient media, such as LB or TY medium, can also be found in the manual by Sambrook et al.

#### Example 1

Preparation of a genomic cosmid gene library from *Corynebacterium glutamicum* ATCC13032

Chromosomal DNA from *Corynebacterium glutamicum* ATCC13032 was isolated as described in Tauch et al. (1995, Plasmid

33, 168-179) and partially cleaved with the restriction enzyme Sau3AI (Amersham Pharmacia, Freiburg, Germany, product description Sau3AI, code no. 27-0913-02). The DNA fragments were dephosphorylated with shrimp alkaline phosphatase (Roche Diagnostics GmbH, Mannheim, Germany, product description SAP, code no. 1758250). The DNA of cosmid vector SuperCos1 (Wahl et al. (1987), Proceedings of the National Academy of Sciences USA 84, 2160-2164), obtained from Stratagene (La Jolla, USA, product description SuperCos1 Cosmid Vector Kit, code no. 251301), was cleaved with the restriction enzyme XbaI (Amersham Pharmacia, Freiburg, Germany, product description XbaI, code no. 27-0948-02) and also dephosphorylated with shrimp alkaline phosphatase.

The cosmid DNA was then cleaved with the restriction enzyme BamHI (Amersham Pharmacia, Freiburg, Germany, product description BamHI, code no. 27-0868-04). The cosmid DNA treated in this way was mixed with the treated ATCC13032 DNA and the mixture was treated with T4 DNA ligase (Amersham Pharmacia, Freiburg, Germany, product description T4 DNA ligase, code no. 27-0870-04). The ligation mixture was then packaged into phages using Gigapack II XL Packing Extract (Stratagene, La Jolla, USA, product description Gigapack II XL Packing Extract, code no. 200217).

For infection of the E. coli strain NM554 (Raleigh et al., 1988, Nucleic Acid Research 16, 1563-1575), the cells were taken up in 10 mM MgSO<sub>4</sub> and mixed with an aliquot of the phage suspension. Infection and titering of the cosmid library were carried out as described in Sambrook et al. (1989, Molecular Cloning: A Laboratory Manual, Cold Spring Harbor), the cells being plated on LB agar (Lennox, 1955, Virology 1, 190) containing 100 mg/l of ampicillin. After incubation overnight at 37°C, recombinant single clones were selected.

## Example 2

### Isolation and sequencing of the pknB gene

The cosmid DNA of a single colony was isolated with the Qiaprep Spin Miniprep Kit (product no. 27106, Qiagen, Hilden, Germany) in accordance with the manufacturer's instructions and partially cleaved with the restriction enzyme Sau3AI (Amersham Pharmacia, Freiburg, Germany, product description Sau3AI, product no. 27-0913-02). The DNA fragments were dephosphorylated with shrimp alkaline phosphatase (Roche Diagnostics GmbH, Mannheim, Germany, product description SAP, product no. 1758250). After separation by gel electrophoresis, the cosmid fragments in the size range from 1500 to 2000 bp were isolated with the QiaExII Gel Extraction Kit (product no. 20021, Qiagen, Hilden, Germany).

The DNA of sequencing vector pZero-1, obtained from Invitrogen (Groningen, The Netherlands, product description Zero Background Cloning Kit, product no. K2500-01), was cleaved with the restriction enzyme BamHI (Amersham Pharmacia, Freiburg, Germany, product description BamHI, product no. 27-0868-04). Ligation of the cosmid fragments into sequencing vector pZero-1 was carried out as described in Sambrook et al. (1989, Molecular Cloning: A Laboratory Manual, Cold Spring Harbor), the DNA mixture being incubated overnight with T4 ligase (Pharmacia Biotech, Freiburg, Germany). This ligation mixture was then introduced into the E. coli strain DH5 $\alpha$ MCR (Grant, 1990, Proceedings of the National Academy of Sciences USA 87, 4645-4649) by electroporation (Tauch et al. 1994, FEMS Microbiol. Letters 123, 343-7) and plated on LB agar (Lennox, 1955, Virology 1, 190) containing 50 mg/l of zeocin.

Plasmid preparation of the recombinant clones was carried out with Biorobot 9600 (product no. 900200, Qiagen, Hilden,

Germany). Sequencing was carried out by the dideoxy chain termination method of Sanger et al. (1977, Proceedings of the National Academy of Sciences USA 74, 5463-5467) with modifications by Zimmermann et al. (1990, Nucleic Acids Research 18, 1067). The "RR dRhodamin Terminator Cycle Sequencing Kit" from PE Applied Biosystems (product no. 403044, Weiterstadt, Germany) was used. Separation by gel electrophoresis and analysis of the sequencing reaction were carried out in a "Rotiphorese NF acrylamide/bisacrylamide" gel (29:1) (product no. A124.1, Roth, Karlsruhe, Germany) with the "ABI Prism 377" sequencer from PE Applied Biosystems (Weiterstadt, Germany).

The raw sequence data obtained were then processed using the Staden programming package (1986, Nucleic Acids Research 14, 217-231), version 97-0. The individual sequences of the pZero-1 derivatives were assembled into a cohesive contig. Computer-assisted coding region analysis was performed with the XNIP program (Staden, 1986, Nucleic Acids Research 14, 217-231).

The nucleotide sequence obtained is shown in SEQ ID No. 1. Analysis of the nucleotide sequence gave an open reading frame of 1884 base pairs, which was called the pknB gene. The pknB gene codes for a protein of 627 amino acids.

**What is claimed is:**

1. An isolated polynucleotide from corynebacteria which contains a polynucleotide sequence coding for the pknB gene and is selected from the group comprising:
  - a) a polynucleotide which is at least 70% identical to a polynucleotide coding for a polypeptide containing the amino acid sequence of SEQ ID No. 2,
  - b) a polynucleotide coding for a polypeptide containing an amino acid sequence which is at least 70% identical to the amino acid sequence of SEQ ID No. 2,
  - c) a polynucleotide which is complementary to the polynucleotides of a) or b), and
  - d) a polynucleotide containing at least 15 consecutive nucleotides of the polynucleotide sequence of a), b) or c),the polypeptide preferably having the activity of protein kinase B.
2. A polynucleotide as claimed in claim 1 which is a preferably recombinant DNA replicatable in corynebacteria.
3. A polynucleotide as claimed in claim 1 which is an RNA.
4. A polynucleotide as claimed in claim 2 which contains the nucleic acid sequence as shown in SEQ ID No. 1.
5. A replicatable DNA as claimed in claim 2 which contains:

- (i) the nucleotide sequence shown in SEQ ID No. 1,  
or
  - (ii) at least one sequence corresponding to sequence  
(i) within the degeneracy of the genetic code,  
or
  - (iii) at least one sequence which hybridizes with the  
sequence complementary to sequence (i) or (ii),  
and optionally
  - (iv) neutral sense mutations in (i).
- 6. A replicatable DNA as claimed in claim 5 wherein the  
hybridization is carried out under a stringency  
corresponding to at most 2x SSC.
  - 7. A polynucleotide sequence as claimed in claim 1 which  
codes for a polypeptide containing the amino acid  
sequence shown in SEQ ID No. 2.
  - 8. Corynebacteria in which the pknB gene is amplified  
and, in particular, overexpressed.
  - 9. A fermentation process for the preparation of L-amino  
acids, especially L-lysine, wherein the following  
steps are carried out:
    - a) fermentation of the corynebacteria producing the  
desired L-amino acid, in which at least the pknB  
gene or nucleotide sequences coding therefor are  
amplified and, in particular, overexpressed,
    - b) enrichment of the L-amino acid in the medium or in  
the cells of the bacteria, and
    - c) isolation of the L-amino acid.
  - 10. The process as claimed in claim 9 wherein bacteria are  
used in which other genes of the biosynthetic pathway

of the desired L-amino acid are additionally amplified.

11. The process as claimed in claim 9 wherein bacteria are used in which the metabolic pathways which reduce the formation of the desired L-amino acid are at least partially switched off.
12. The process as claimed in claim 9 wherein a strain transformed with a plasmid vector is used and the plasmid vector carries the nucleotide sequence coding for the *pknB* gene.
13. The process as claimed in claim 9 wherein the expression of the polynucleotide(s) coding for the *pknB* gene is amplified and, in particular, overexpressed.
14. The process as claimed in claim 9 wherein the catalytic properties of the polypeptide (enzyme protein) for which the *pknB* polynucleotide codes are enhanced.
15. The process as claimed in claim 9 wherein, for the production of L-amino acids, coryneform microorganisms are fermented in which one or more endogeny genes selected from the following group are simultaneously amplified or overexpressed:
  - 15.1 the *dapA* gene coding for dihydrodipicolinate synthase,
  - 15.2 the *gap* gene coding for glyceraldehyde 3-phosphate dehydrogenase,
  - 15.3 the *tpi* gene coding for triose phosphate isomerase,
  - 15.4 the *pgk* gene coding for 3-phosphoglycerate kinase,

- 15.5 the zwf gene coding for glucose-6-phosphate dehydrogenase,
  - 15.6 the pyc gene coding for pyruvate carboxylase,
  - 15.7 the lysC gene coding for a feedback-resistant aspartate kinase,
  - 15.8 the lysE gene coding for lysine export,
  - 15.9 the hom gene coding for homoserine dehydrogenase,
  - 15.10 the ilvA gene coding for threonine dehydratase or the ilvA(Fbr) allele coding for a feedback-resistant threonine dehydratase,
  - 15.11 the ilvBN gene coding for acetohydroxy acid synthase,
  - 15.12 the ilvD gene coding for dihydroxy acid dehydratase, or
  - 15.13 the zwal gene coding for the Zwal protein.
16. The process as claimed in claim 9 wherein, for the production of L-amino acids, coryneform microorganisms are fermented in which one or more genes selected from the following group are simultaneously attenuated:
- 16.1 the pck gene coding for phosphoenol pyruvate carboxykinase,
  - 16.2 the pgi gene coding for glucose-6-phosphate isomerase,
  - 16.3 the poxB gene coding for pyruvate oxidase, or
  - 16.4 the zwa2 gene coding for the Zwa2 protein.

17. Corynebacteria which contain a vector carrying a polynucleotide as claimed in claim 1.
18. The process as claimed in one or more of claims 9-16, wherein microorganisms of the species *Corynebacterium glutamicum* are used.
19. The process as claimed in claim 18, wherein the *Corynebacterium glutamicum* strain DSM 13994 is employed.
20. A method of detecting RNA, cDNA and DNA in order to isolate nucleic acids, or polynucleotides or genes, which code for protein kinase B or have a high degree of similarity to the sequence of the pknB gene, wherein the polynucleotide containing the polynucleotide sequences as claimed in claims 1, 2, 3 or 4 is used as hybridization probes.
21. The method as claimed in claim 20 wherein arrays, micro-arrays or DNA chips are used.
22. A DNA originating from corynebacteria and coding for protein kinase B, wherein the corresponding amino acid sequences between positions 581 and 587 in SEQ ID No. 2 are modified by amino acid exchange.
23. A DNA originating from corynebacteria and coding for protein kinase B, wherein the corresponding amino acid sequences contain any other proteogenic amino acid except L-proline in position 584 in SEQ ID No. 2.
24. A DNA originating from corynebacteria and coding for protein kinase B, wherein the corresponding amino acid sequences contain L-serine or L-threonine in position 584 in SEQ ID No. 2.
25. A DNA originating from corynebacteria and coding for protein kinase B, wherein the corresponding amino acid

sequence contains L-serine in position 584, shown in SEQ ID No. 4.

26. A DNA as claimed in claim 25 which contains the nucleobase thymine in position 2343, shown in SEQ ID No. 3.
27. Corynebacteria which contain a DNA as claimed in claim 22, 23, 24, 25 or 26.
28. Corynebacterium glutamicum DM1547 deposited as DSM13994 in the Deutsche Sammlung für Mikroorganismen und Zellkulturen (German Collection of Microorganisms and Cell Cultures), Brunswick, Germany.

## SEQUENCE LISTING

&lt;110&gt; Degussa AG

&lt;120&gt; Nucleotide sequences coding for the pknB gene

&lt;130&gt; 000505 BT

&lt;140&gt;

&lt;141&gt;

&lt;160&gt; 4

&lt;170&gt; PatentIn Ver. 2.1

&lt;210&gt; 1

&lt;211&gt; 2875

&lt;212&gt; DNA

&lt;213&gt; Corynebacterium glutamicum

&lt;220&gt;

&lt;221&gt; CDS

&lt;222&gt; (594) .. (2474)

&lt;223&gt; pknB gene

&lt;400&gt; 1

```

cgcacgagcg cgatgatggc gcaggcggag gcgccgtcgc caagcgaatc aacggcgatg 60
ctgggcaggg tggcccggcc tgcaacaatc acccaagaag cggccccgaa acgcggttcc 120
ggcattggca ttggtctgtt catcgcagct ttgcttgccg tgattattgg cgcggtgatc 180
tatgcgggca ccaccggaat tttgttcaac gacactccgg aagaaaccac cacacctgaa 240
accattacgg aacatacac cccaaccgtg gaggaacca cctctcagtg ggtaccgcca 300
acgcctccaa cacggtcaac attcaccgaa cctgaaacaa cttcacaccg tccgacgaca 360
agtgaagaga gcacatccga ggaaccaacc acggaagctc caacaagtag ccgaactgtg 420
cctcaaatcc ctacctctac acctaggacg agtgctagcg ttccagttga gactaatgca 480
ccggctgatg atttaatcga cgccgtaaat ggcctattgg atgtaggagg agcgcagtga 540
ccttcgtgat cgctgatcgc tatgaactgg atgccgtcat cggctccggt ggc atg 596
                                     Met
                                     1

agc gag gtg ttc gcg gcc acc gac acg ctc att ggt cgg gag gtc gcg 644
Ser Glu Val Phe Ala Ala Thr Asp Thr Leu Ile Gly Arg Glu Val Ala
                    5                      10                      15

gta aag atg ctg cgc atc gac ctt gcg aaa gat ccc aat ttc cga gaa 692
Val Lys Met Leu Arg Ile Asp Leu Ala Lys Asp Pro Asn Phe Arg Glu
                    20                      25                      30

cgc ttc cgc agg gaa gcc caa aac tcc gga agg ttg agc cac tct tcg 740
Arg Phe Arg Arg Glu Ala Gln Asn Ser Gly Arg Leu Ser His Ser Ser
                    35                      40                      45

```

atc gtc gct gtt ttt gac acc ggc gaa gta gac aaa gac ggc acc tct	788
Ile Val Ala Val Phe Asp Thr Gly Glu Val Asp Lys Asp Gly Thr Ser	
50 55 60 65	
gtt ccc tac att gtg atg gaa cgc gtg cag ggt cga aac ctg cgc gaa	836
Val Pro Tyr Ile Val Met Glu Arg Val Gln Gly Arg Asn Leu Arg Glu	
70 75 80	
gtt gtc acc gaa gac ggc gta ttc acc cca gtt gag gca gcc aac atc	884
Val Val Thr Glu Asp Gly Val Phe Thr Pro Val Glu Ala Ala Asn Ile	
85 90 95	
ctc atc cct gtg tgt gaa gcg ctg cag gca tcc cat gac gcc ggc att	932
Leu Ile Pro Val Cys Glu Ala Leu Gln Ala Ser His Asp Ala Gly Ile	
100 105 110	
att cac cgc gat gtg aaa ccc gcc aac atc atg atc acc aac acc ggt	980
Ile His Arg Asp Val Lys Pro Ala Asn Ile Met Ile Thr Asn Thr Gly	
115 120 125	
ggc gtg aaa gtc atg gac ttc ggc atc gcc cgc gcg gtc aac gat tcc	1028
Gly Val Lys Val Met Asp Phe Gly Ile Ala Arg Ala Val Asn Asp Ser	
130 135 140 145	
acc tcc gcc atg act caa acc tcc gca gtc atc ggc acc gcc cag tac	1076
Thr Ser Ala Met Thr Gln Thr Ser Ala Val Ile Gly Thr Ala Gln Tyr	
150 155 160	
ctc tcc cct gag cag gcc cgc ggc aaa ccc gcc gat gcg cgt tcc gat	1124
Leu Ser Pro Glu Gln Ala Arg Gly Lys Pro Ala Asp Ala Arg Ser Asp	
165 170 175	
att tac gcc acc ggc tgc gtc atg tac gaa tta gtc acc ggt aag cca	1172
Ile Tyr Ala Thr Gly Cys Val Met Tyr Glu Leu Val Thr Gly Lys Pro	
180 185 190	
cct ttt gaa ggc gag tcc cct ttc gcc gtg gcc tac caa cac gtc cag	1220
Pro Phe Glu Gly Glu Ser Pro Phe Ala Val Ala Tyr Gln His Val Gln	
195 200 205	
gaa gac ccc acc cct cct tcg gat ttc atc gcg gac ctc acc ccg acc	1268
Glu Asp Pro Thr Pro Pro Ser Asp Phe Ile Ala Asp Leu Thr Pro Thr	
210 215 220 225	
tct gct gtc aac gtg gat gcc gtg gta ctc acc gcc atg gca aaa cac	1316
Ser Ala Val Asn Val Asp Ala Val Val Leu Thr Ala Met Ala Lys His	
230 235 240	
ccc gcc gac cgc tac caa aca gcc tcc gaa atg gcc gct gac ctg ggc	1364
Pro Ala Asp Arg Tyr Gln Thr Ala Ser Glu Met Ala Ala Asp Leu Gly	
245 250 255	
cgg cta tcc cgc aat gca gtc tcc cat gcc gca cgc gcg cat gta gaa	1412
Arg Leu Ser Arg Asn Ala Val Ser His Ala Ala Arg Ala His Val Glu	
260 265 270	
aca gaa gaa acc cca gaa gag ccc gaa act cgc ttc tcg acg cgc acc	1460
Thr Glu Glu Thr Pro Glu Glu Pro Glu Thr Arg Phe Ser Thr Arg Thr	
275 280 285	

tcc acc caa gtg gcc ccc gcc gca ggc gtg gct gcg gcc agt acg ggg Ser Thr Gln Val Ala Pro Ala Ala Gly Val Ala Ala Ala Ser Thr Gly 290 295 300 305	1508
tca ggg tct tct tcg cgt aaa cgt gga tcc aga ggc ctc acc gcc ctg Ser Gly Ser Ser Arg Lys Arg Gly Ser Arg Gly Leu Thr Ala Leu 310 315 320	1556
gcc atc gtg tta tcc cta ggt gtc gtc ggc gtt gcc ggt gcc ttc acc Ala Ile Val Leu Ser Leu Gly Val Val Gly Val Ala Gly Ala Phe Thr 325 330 335	1604
tac gac tac ttt gcc aac agc tcc tcc act gca acc agc gcg atc ccc Tyr Asp Tyr Phe Ala Asn Ser Ser Ser Thr Ala Thr Ser Ala Ile Pro 340 345 350	1652
aat gtg gaa ggc ctc ccg cag caa gaa gct ctc aca gaa ctt caa gca Asn Val Glu Gly Leu Pro Gln Gln Glu Ala Leu Thr Glu Leu Gln Ala 355 360 365	1700
gca gga ttt gtt gtc aac atc gtc gaa gaa gcc agc gcc gac gtc gcc Ala Gly Phe Val Val Asn Ile Val Glu Glu Ala Ser Ala Asp Val Ala 370 375 380 385	1748
gaa ggc ctc gtc atc cga gca aac cca agc gtt gga tcc gaa atc cgc Glu Gly Leu Val Ile Arg Ala Asn Pro Ser Val Gly Ser Glu Ile Arg 390 395 400	1796
caa ggg gcc acc gtc acc atc acc gtg tcc acc ggc cga gaa atg atc Gln Gly Ala Thr Val Thr Ile Thr Val Ser Thr Gly Arg Glu Met Ile 405 410 415	1844
aac atc cca gac gtc tcc ggc atg aca ctt gag gac gcc gcc cgc gcc Asn Ile Pro Asp Val Ser Gly Met Thr Leu Glu Asp Ala Ala Arg Ala 420 425 430	1892
ctc gaa gac gtt ggt ctc ata ctc aac caa aac gtt cgg gaa gaa acc Leu Glu Asp Val Gly Leu Ile Leu Asn Gln Asn Val Arg Glu Glu Thr 435 440 445	1940
tcc gac gac gtc gaa tct ggc ctc gtc atc gac caa aac ccc gaa gcc Ser Asp Asp Val Glu Ser Gly Leu Val Ile Asp Gln Asn Pro Glu Ala 450 455 460 465	1988
ggc caa gaa gta gtc gtg ggt tcc tct gta tct cta acc atg tct tca Gly Gln Glu Val Val Val Gly Ser Ser Val Ser Leu Thr Met Ser Ser 470 475 480	2036
ggc acc gag agc atc cga gtg ccc aac ctc acc ggc atg aac tgg tca Gly Thr Glu Ser Ile Arg Val Pro Asn Leu Thr Gly Met Asn Trp Ser 485 490 495	2084
caa gca gaa caa aac ctc atc tcc atg ggc ttt aac ccc aca gct tcc Gln Ala Glu Gln Asn Leu Ile Ser Met Gly Phe Asn Pro Thr Ala Ser 500 505 510	2132
tac tta gac agc agc gaa cca gaa ggc gaa gtc ctc tca gtt tcc agc Tyr Leu Asp Ser Ser Glu Pro Glu Gly Glu Val Leu Ser Val Ser Ser 515 520 525	2180

caa gga act gaa cta ccc aag ggt tca tcc atc aca gtg gaa gtc tcc 2228  
 Gln Gly Thr Glu Leu Pro Lys Gly Ser Ser Ile Thr Val Glu Val Ser  
 530 535 540 545

aac ggc atg ctc atc caa gcc ccc gat ctc gcc cgc atg tcc acc gaa 2276  
 Asn Gly Met Leu Ile Gln Ala Pro Asp Leu Ala Arg Met Ser Thr Glu  
 550 555 560

cag gcc atc agt gcc ctc cgc gct gct ggc tgg acc gcc cca gat caa 2324  
 Gln Ala Ile Ser Ala Leu Arg Ala Gly Trp Thr Ala Pro Asp Gln  
 565 570 575

tcc ctg atc gtc ggc gac ccc atc cac acc gca gcc ctc gtg gat caa 2372  
 Ser Leu Ile Val Gly Asp Pro Ile His Thr Ala Ala Leu Val Asp Gln  
 580 585 590

aac aaa atc gga ttc caa tcc cca acc cct gca acc ctc ttc cgc aaa 2420  
 Asn Lys Ile Gly Phe Gln Ser Pro Thr Pro Ala Thr Leu Phe Arg Lys  
 595 600 605

gac gcc caa gtg caa gtg cga ctc ttc gaa ttc gat ctc gct gca ctc 2468  
 Asp Ala Gln Val Gln Val Arg Leu Phe Glu Phe Asp Leu Ala Ala Leu  
 610 615 620 625

gtg caa tagccaacaa ggaaaccgtc aaggtagctg gcccggaac tgatacgtta 2524  
 Val Gln

agctcaaaca agataagtac cagttgctgg ggtttttcca agacaataaa ttatgaaggt 2584

gtgaacaatg ccaaaggcaa gagtaactaa aaacgagacc gcaccggttt caagcaaccc 2644

aagcgcaaac cgacccccgg ttaagatcaa ttccgcccga accccaatgt ggtacaaggt 2704

catcatgttt gccttcatga tcgtcggcct agcctggttg atcattaact acctcgtggg 2764

cccacagatc ccattcatgg ctgatcttgg tgcattggaac tatggcatcg gcttcggtct 2824

gatgatcatc ggcctactca tgaccatggg ttggcggttaa tccttcaaaa a 2875

<210> 2  
 <211> 627  
 <212> PRT  
 <213> Corynebacterium glutamicum

<400> 2  
 Met Ser Glu Val Phe Ala Ala Thr Asp Thr Leu Ile Gly Arg Glu Val  
 1 5 10 15

Ala Val Lys Met Leu Arg Ile Asp Leu Ala Lys Asp Pro Asn Phe Arg  
 20 25 30

Glu Arg Phe Arg Arg Glu Ala Gln Asn Ser Gly Arg Leu Ser His Ser  
 35 40 45

Ser Ile Val Ala Val Phe Asp Thr Gly Glu Val Asp Lys Asp Gly Thr  
 50 55 60

Ser Val Pro Tyr Ile Val Met Glu Arg Val Gln Gly Arg Asn Leu Arg  
 65 70 75 80  
 Glu Val Val Thr Glu Asp Gly Val Phe Thr Pro Val Glu Ala Ala Asn  
 85 90 95  
 Ile Leu Ile Pro Val Cys Glu Ala Leu Gln Ala Ser His Asp Ala Gly  
 100 105 110  
 Ile Ile His Arg Asp Val Lys Pro Ala Asn Ile Met Ile Thr Asn Thr  
 115 120 125  
 Gly Gly Val Lys Val Met Asp Phe Gly Ile Ala Arg Ala Val Asn Asp  
 130 135 140  
 Ser Thr Ser Ala Met Thr Gln Thr Ser Ala Val Ile Gly Thr Ala Gln  
 145 150 155 160  
 Tyr Leu Ser Pro Glu Gln Ala Arg Gly Lys Pro Ala Asp Ala Arg Ser  
 165 170 175  
 Asp Ile Tyr Ala Thr Gly Cys Val Met Tyr Glu Leu Val Thr Gly Lys  
 180 185 190  
 Pro Pro Phe Glu Gly Glu Ser Pro Phe Ala Val Ala Tyr Gln His Val  
 195 200 205  
 Gln Glu Asp Pro Thr Pro Pro Ser Asp Phe Ile Ala Asp Leu Thr Pro  
 210 215 220  
 Thr Ser Ala Val Asn Val Asp Ala Val Val Leu Thr Ala Met Ala Lys  
 225 230 235 240  
 His Pro Ala Asp Arg Tyr Gln Thr Ala Ser Glu Met Ala Ala Asp Leu  
 245 250 255  
 Gly Arg Leu Ser Arg Asn Ala Val Ser His Ala Ala Arg Ala His Val  
 260 265 270  
 Glu Thr Glu Glu Thr Pro Glu Glu Pro Glu Thr Arg Phe Ser Thr Arg  
 275 280 285  
 Thr Ser Thr Gln Val Ala Pro Ala Ala Gly Val Ala Ala Ala Ser Thr  
 290 295 300  
 Gly Ser Gly Ser Ser Ser Arg Lys Arg Gly Ser Arg Gly Leu Thr Ala  
 305 310 315 320  
 Leu Ala Ile Val Leu Ser Leu Gly Val Val Gly Val Ala Gly Ala Phe  
 325 330 335  
 Thr Tyr Asp Tyr Phe Ala Asn Ser Ser Ser Thr Ala Thr Ser Ala Ile  
 340 345 350  
 Pro Asn Val Glu Gly Leu Pro Gln Gln Glu Ala Leu Thr Glu Leu Gln  
 355 360 365  
 Ala Ala Gly Phe Val Val Asn Ile Val Glu Glu Ala Ser Ala Asp Val  
 370 375 380

Ala Glu Gly Leu Val Ile Arg Ala Asn Pro Ser Val Gly Ser Glu Ile  
 385 390 395 400  
 Arg Gln Gly Ala Thr Val Thr Ile Thr Val Ser Thr Gly Arg Glu Met  
 405 410 415  
 Ile Asn Ile Pro Asp Val Ser Gly Met Thr Leu Glu Asp Ala Ala Arg  
 420 425 430  
 Ala Leu Glu Asp Val Gly Leu Ile Leu Asn Gln Asn Val Arg Glu Glu  
 435 440 445  
 Thr Ser Asp Asp Val Glu Ser Gly Leu Val Ile Asp Gln Asn Pro Glu  
 450 455 460  
 Ala Gly Gln Glu Val Val Val Gly Ser Ser Val Ser Leu Thr Met Ser  
 465 470 475 480  
 Ser Gly Thr Glu Ser Ile Arg Val Pro Asn Leu Thr Gly Met Asn Trp  
 485 490 495  
 Ser Gln Ala Glu Gln Asn Leu Ile Ser Met Gly Phe Asn Pro Thr Ala  
 500 505 510  
 Ser Tyr Leu Asp Ser Ser Glu Pro Glu Gly Glu Val Leu Ser Val Ser  
 515 520 525  
 Ser Gln Gly Thr Glu Leu Pro Lys Gly Ser Ser Ile Thr Val Glu Val  
 530 535 540  
 Ser Asn Gly Met Leu Ile Gln Ala Pro Asp Leu Ala Arg Met Ser Thr  
 545 550 555 560  
 Glu Gln Ala Ile Ser Ala Leu Arg Ala Ala Gly Trp Thr Ala Pro Asp  
 565 570 575  
 Gln Ser Leu Ile Val Gly Asp Pro Ile His Thr Ala Ala Leu Val Asp  
 580 585 590  
 Gln Asn Lys Ile Gly Phe Gln Ser Pro Thr Pro Ala Thr Leu Phe Arg  
 595 600 605  
 Lys Asp Ala Gln Val Gln Val Arg Leu Phe Glu Phe Asp Leu Ala Ala  
 610 615 620  
 Leu Val Gln  
 625

<210> 3  
 <211> 2875  
 <212> DNA  
 <213> Corynebacterium glutamicum

<220>  
 <221> CDS  
 <222> (594)..(2474)

<223> pknB-1547 allele

<220>

<221> misc feature

<222> (2343)

<223> C-T Transition

<400> 3

```

cgcacgagcg cgatgatggc gcaggcggag gcgccgtcgc caagcgaatc aacggcgatg 60
ctgggcaggg tggcccggcc tgcaacaatc acccaagaag cggccccgaa acgcgggttcc 120
ggcattggca ttggtctggt catcgagct ttgcttgccg tgattattgg cgcgggtgatc 180
tatgcgggca ccaccggaat tttgttcaac gacactccgg aagaaaccac cacacctgaa 240
accattacgg aacatacac cccaaccgtg gaggaacca cctctcagtg ggtaccgcca 300
acgcctccaa cacgggtcaac attcaccgaa cctgaaacaa cttcacaccg tccgacgaca 360
agtgaagaga gcacatccga ggaaccaacc acggaagctc caacaagtag ccgaactgtg 420
cctcaaatcc ctacctctac acctaggacg agtgctagcg ttccagttga gactaatgca 480
ccggctgatg atttaatcga cgccgtaaat ggcctattgg atgtaggagg agcgcagtga 540
ccttcgtgat cgctgatcgc tatgaactgg atgccgtcat cggctccggt ggc atg 596
                                     Met
                                     1

agc gag gtg ttc gcg gcc acc gac acg ctc att ggt cgg gag gtc gcg 644
Ser Glu Val Phe Ala Ala Thr Asp Thr Leu Ile Gly Arg Glu Val Ala
                                     5
                                     10
                                     15

gta aag atg ctg cgc atc gac ctt gcg aaa gat ccc aat ttc cga gaa 692
Val Lys Met Leu Arg Ile Asp Leu Ala Lys Asp Pro Asn Phe Arg Glu
                                     20
                                     25
                                     30

cgc ttc cgc agg gaa gcc caa aac tcc gga agg ttg agc cac tct tcg 740
Arg Phe Arg Arg Glu Ala Gln Asn Ser Gly Arg Leu Ser His Ser Ser
                                     35
                                     40
                                     45

atc gtc gct gtt ttt gac acc ggc gaa gta gac aaa gac ggc acc tct 788
Ile Val Ala Val Phe Asp Thr Gly Glu Val Asp Lys Asp Gly Thr Ser
                                     50
                                     55
                                     60
                                     65

gtt ccc tac att gtg atg gaa cgc gtg cag ggt cga aac ctg cgc gaa 836
Val Pro Tyr Ile Val Met Glu Arg Val Gln Gly Arg Asn Leu Arg Glu
                                     70
                                     75
                                     80

gtt gtc acc gaa gac ggc gta ttc acc cca gtt gag gca gcc aac atc 884
Val Val Thr Glu Asp Gly Val Phe Thr Pro Val Glu Ala Ala Asn Ile
                                     85
                                     90
                                     95

ctc atc cct gtg tgt gaa gcg ctg cag gca tcc cat gac gcc ggc att 932
Leu Ile Pro Val Cys Glu Ala Leu Gln Ala Ser His Asp Ala Gly Ile
                                     100
                                     105
                                     110

```

att cac cgc gat gtg aaa ccc gcc aac atc atg atc acc aac acc ggt Ile His Arg Asp Val Lys Pro Ala Asn Ile Met Ile Thr Asn Thr Gly 115 120 125	980
ggc gtg aaa gtc atg gac ttc ggc atc gcc cgc gcg gtc aac gat tcc Gly Val Lys Val Met Asp Phe Gly Ile Ala Arg Ala Val Asn Asp Ser 130 135 140 145	1028
acc tcc gcc atg act caa acc tcc gca gtc atc ggc acc gcc cag tac Thr Ser Ala Met Thr Gln Thr Ser Ala Val Ile Gly Thr Ala Gln Tyr 150 155 160	1076
ctc tcc cct gag cag gcc cgc ggc aaa ccc gcc gat gcg cgt tcc gat Leu Ser Pro Glu Gln Ala Arg Gly Lys Pro Ala Asp Ala Arg Ser Asp 165 170 175	1124
att tac gcc acc ggc tgc gtc atg tac gaa tta gtc acc ggt aag cca Ile Tyr Ala Thr Gly Cys Val Met Tyr Glu Leu Val Thr Gly Lys Pro 180 185 190	1172
cct ttt gaa ggc gag tcc cct ttc gcc gtg gcc tac caa cac gtc cag Pro Phe Glu Gly Glu Ser Pro Phe Ala Val Ala Tyr Gln His Val Gln 195 200 205	1220
gaa gac ccc acc cct cct tcg gat ttc atc gcg gac ctc acc ccg acc Glu Asp Pro Thr Pro Pro Ser Asp Phe Ile Ala Asp Leu Thr Pro Thr 210 215 220 225	1268
tct gct gtc aac gtg gat gcc gtg gta ctc acc gcc atg gca aaa cac Ser Ala Val Asn Asp Ala Val Val Leu Thr Ala Met Ala Lys His 230 235 240	1316
ccc gcc gac cgc tac caa aca gcc tcc gaa atg gcc gct gac ctg ggc Pro Ala Asp Arg Tyr Gln Thr Ala Ser Glu Met Ala Ala Asp Leu Gly 245 250 255	1364
cgg cta tcc cgc aat gca gtc tcc cat gcc gca cgc gcg cat gta gaa Arg Leu Ser Arg Asn Ala Val Ser His Ala Ala Arg Ala His Val Glu 260 265 270	1412
aca gaa gaa acc cca gaa gag ccc gaa act cgc ttc tcg acg cgc acc Thr Glu Glu Thr Pro Glu Glu Pro Glu Thr Arg Phe Ser Thr Arg Thr 275 280 285	1460
tcc acc caa gtg gcc ccc gcc gca ggc gtg gct gcg gcc agt acg ggg Ser Thr Gln Val Ala Pro Ala Ala Gly Val Ala Ala Ser Thr Gly 290 295 300 305	1508
tca ggg tct tct tcg cgt aaa cgt gga tcc aga ggc ctc acc gcc ctg Ser Gly Ser Ser Ser Arg Lys Arg Gly Ser Arg Gly Leu Thr Ala Leu 310 315 320	1556
gcc atc gtg tta tcc cta ggt gtc gtc ggc gtt gcc ggt gcc ttc acc Ala Ile Val Leu Ser Leu Gly Val Val Gly Val Ala Gly Ala Phe Thr 325 330 335	1604
tac gac tac ttt gcc aac agc tcc tcc act gca acc agc gcg atc ccc Tyr Asp Tyr Phe Ala Asn Ser Ser Ser Thr Ala Thr Ser Ala Ile Pro 340 345 350	1652

aat gtg gaa ggc ctc ccg cag caa gaa gct ctc aca gaa ctt caa gca	1700
Asn Val Glu Gly Leu Pro Gln Gln Glu Ala Leu Thr Glu Leu Gln Ala	
355 360 365	
gca gga ttt gtt gtc aac atc gtc gaa gaa gcc agc gcc gac gtc gcc	1748
Ala Gly Phe Val Val Asn Ile Val Glu Glu Ala Ser Ala Asp Val Ala	
370 375 380 385	
gaa ggc ctc gtc atc cga gca aac cca agc gtt gga tcc gaa atc cgc	1796
Glu Gly Leu Val Ile Arg Ala Asn Pro Ser Val Gly Ser Glu Ile Arg	
390 395 400	
caa ggg gcc acc gtc acc atc acc gtg tcc acc ggc cga gaa atg atc	1844
Gln Gly Ala Thr Val Thr Ile Thr Val Ser Thr Gly Arg Glu Met Ile	
405 410 415	
aac atc cca gac gtc tcc ggc atg aca ctt gag gac gcc gcc cgc gcc	1892
Asn Ile Pro Asp Val Ser Gly Met Thr Leu Glu Asp Ala Ala Arg Ala	
420 425 430	
ctc gaa gac gtt ggt ctc ata ctc aac caa aac gtt cgg gaa gaa acc	1940
Leu Glu Asp Val Gly Leu Ile Leu Asn Gln Asn Val Arg Glu Glu Thr	
435 440 445	
tcc gac gac gtc gaa tct ggc ctc gtc atc gac caa aac ccc gaa gcc	1988
Ser Asp Asp Val Glu Ser Gly Leu Val Ile Asp Gln Asn Pro Glu Ala	
450 455 460 465	
ggc caa gaa gta gtc gtg ggt tcc tct gta tct cta acc atg tct tca	2036
Gly Gln Glu Val Val Val Gly Ser Ser Val Ser Leu Thr Met Ser Ser	
470 475 480	
ggc acc gag agc atc cga gtg ccc aac ctc acc ggc atg aac tgg tca	2084
Gly Thr Glu Ser Ile Arg Val Pro Asn Leu Thr Gly Met Asn Trp Ser	
485 490 495	
caa gca gaa caa aac ctc atc tcc atg ggc ttt aac ccc aca gct tcc	2132
Gln Ala Glu Gln Asn Leu Ile Ser Met Gly Phe Asn Pro Thr Ala Ser	
500 505 510	
tac tta gac agc agc gaa cca gaa ggc gaa gtc ctc tca gtt tcc agc	2180
Tyr Leu Asp Ser Ser Glu Pro Glu Gly Glu Val Leu Ser Val Ser Ser	
515 520 525	
caa gga act gaa cta ccc aag ggt tca tcc atc aca gtg gaa gtc tcc	2228
Gln Gly Thr Glu Leu Pro Lys Gly Ser Ser Ile Thr Val Glu Val Ser	
530 535 540 545	
aac ggc atg ctc atc caa gcc ccc gat ctc gcc cgc atg tcc acc gaa	2276
Asn Gly Met Leu Ile Gln Ala Pro Asp Leu Ala Arg Met Ser Thr Glu	
550 555 560	
cag gcc atc agt gcc ctc cgc gct gct ggc tgg acc gcc cca gat caa	2324
Gln Ala Ile Ser Ala Leu Arg Ala Ala Gly Trp Thr Ala Pro Asp Gln	
565 570 575	
tcc ctg atc gtc ggc gac tcc atc cac acc gca gcc ctc gtg gat caa	2372
Ser Leu Ile Val Gly Asp Ser Ile His Thr Ala Ala Leu Val Asp Gln	
580 585 590	

aac aaa atc gga ttc caa tcc cca acc cct gca acc ctc ttc cgc aaa 2420  
 Asn Lys Ile Gly Phe Gln Ser Pro Thr Pro Ala Thr Leu Phe Arg Lys  
 595 600 605

gac gcc caa gtg caa gtg cga ctc ttc gaa ttc gat ctc gct gca ctc 2468  
 Asp Ala Gln Val Gln Val Arg Leu Phe Glu Phe Asp Leu Ala Ala Leu  
 610 615 620 625

gtg caa tagccaacaa-ggaaaccgtc-aaggtagctg gcccggaac-tgatacgtta 2524  
 Val Gln

agctcaaaca agataagtac cagttgctgg gggtttttcca agacaataaa ttatgaaggt 2584

gtgaacaatg ccaaaggcaa gagtaactaa aaacgagacc gcaccggttt caagcaaccc 2644

aagcgcaaac cgcacccccg ttaagatcaa ttccgccgga accccaatgt ggtacaaggt 2704

catcatgttt gccttcatga tcgtcggcct agcctgggtg atcattaact acctcgtggg 2764

cccacagatc ccattcatgg ctgatcttgg tgcattggaac tatggcatcg gcttcggtct 2824

gatgatcatc ggctactca tgaccatggg ttggcggtta tccttcaaaa a 2875

<210> 4

<211> 627

<212> PRT

<213> Corynebacterium glutamicum

<400> 4

Met Ser Glu Val Phe Ala Ala Thr Asp Thr Leu Ile Gly Arg Glu Val 15  
 1 5 10

Ala Val Lys Met Leu Arg Ile Asp Leu Ala Lys Asp Pro Asn Phe Arg 30  
 20 25 30

Glu Arg Phe Arg Arg Glu Ala Gln Asn Ser Gly Arg Leu Ser His Ser 45  
 35 40 45

Ser Ile Val Ala Val Phe Asp Thr Gly Glu Val Asp Lys Asp Gly Thr 60  
 50 55 60

Ser Val Pro Tyr Ile Val Met Glu Arg Val Gln Gly Arg Asn Leu Arg 80  
 65 70 75 80

Glu Val Val Thr Glu Asp Gly Val Phe Thr Pro Val Glu Ala Ala Asn 95  
 85 90 95

Ile Leu Ile Pro Val Cys Glu Ala Leu Gln Ala Ser His Asp Ala Gly 110  
 100 105 110

Ile Ile His Arg Asp Val Lys Pro Ala Asn Ile Met Ile Thr Asn Thr 125  
 115 120 125

Gly Gly Val Lys Val Met Asp Phe Gly Ile Ala Arg Ala Val Asn Asp 140  
 130 135 140

Ser Thr Ser Ala Met Thr Gln Thr Ser Ala Val Ile Gly Thr Ala Gln  
 145 150 155 160  
 Tyr Leu Ser Pro Glu Gln Ala Arg Gly Lys Pro Ala Asp Ala Arg Ser  
 165 170 175  
 Asp Ile Tyr Ala Thr Gly Cys Val Met Tyr Glu Leu Val Thr Gly Lys  
 180 185 190  
 Pro Pro Phe Glu Gly Glu Ser Pro Phe Ala Val Ala Tyr Gln His Val  
 195 200 205  
 Gln Glu Asp Pro Thr Pro Pro Ser Asp Phe Ile Ala Asp Leu Thr Pro  
 210 215 220  
 Thr Ser Ala Val Asn Val Asp Ala Val Val Leu Thr Ala Met Ala Lys  
 225 230 235 240  
 His Pro Ala Asp Arg Tyr Gln Thr Ala Ser Glu Met Ala Ala Asp Leu  
 245 250 255  
 Gly Arg Leu Ser Arg Asn Ala Val Ser His Ala Ala Arg Ala His Val  
 260 265 270  
 Glu Thr Glu Glu Thr Pro Glu Glu Pro Glu Thr Arg Phe Ser Thr Arg  
 275 280 285  
 Thr Ser Thr Gln Val Ala Pro Ala Ala Gly Val Ala Ala Ala Ser Thr  
 290 295 300  
 Gly Ser Gly Ser Ser Ser Arg Lys Arg Gly Ser Arg Gly Leu Thr Ala  
 305 310 315 320  
 Leu Ala Ile Val Leu Ser Leu Gly Val Val Gly Val Ala Gly Ala Phe  
 325 330 335  
 Thr Tyr Asp Tyr Phe Ala Asn Ser Ser Ser Thr Ala Thr Ser Ala Ile  
 340 345 350  
 Pro Asn Val Glu Gly Leu Pro Gln Gln Glu Ala Leu Thr Glu Leu Gln  
 355 360 365  
 Ala Ala Gly Phe Val Val Asn Ile Val Glu Glu Ala Ser Ala Asp Val  
 370 375 380  
 Ala Glu Gly Leu Val Ile Arg Ala Asn Pro Ser Val Gly Ser Glu Ile  
 385 390 395 400  
 Arg Gln Gly Ala Thr Val Thr Ile Thr Val Ser Thr Gly Arg Glu Met  
 405 410 415  
 Ile Asn Ile Pro Asp Val Ser Gly Met Thr Leu Glu Asp Ala Ala Arg  
 420 425 430  
 Ala Leu Glu Asp Val Gly Leu Ile Leu Asn Gln Asn Val Arg Glu Glu  
 435 440 445  
 Thr Ser Asp Asp Val Glu Ser Gly Leu Val Ile Asp Gln Asn Pro Glu  
 450 455 460

Ala Gly Gln Glu Val Val Val Gly Ser Ser Val Ser Leu Thr Met Ser  
465 470 475 480

Ser Gly Thr Glu Ser Ile Arg Val Pro Asn Leu Thr Gly Met Asn Trp  
485 490 495

Ser Gln Ala Glu Gln Asn Leu Ile Ser Met Gly Phe Asn Pro Thr Ala  
500 505 510

Ser Tyr Leu Asp Ser Ser Glu Pro Glu Gly Glu Val Leu Ser Val Ser  
515 520 525

Ser Gln Gly Thr Glu Leu Pro Lys Gly Ser Ser Ile Thr Val Glu Val  
530 535 540

Ser Asn Gly Met Leu Ile Gln Ala Pro Asp Leu Ala Arg Met Ser Thr  
545 550 555 560

Glu Gln Ala Ile Ser Ala Leu Arg Ala Ala Gly Trp Thr Ala Pro Asp  
565 570 575

Gln Ser Leu Ile Val Gly Asp Ser Ile His Thr Ala Ala Leu Val Asp  
580 585 590

Gln Asn Lys Ile Gly Phe Gln Ser Pro Thr Pro Ala Thr Leu Phe Arg  
595 600 605


Lys Asp Ala Gln Val Gln Val Arg Leu Phe Glu Phe Asp Leu Ala Ala  
610 615 620

Leu Val Gln  
625

0-1	Form - PCT/RO/134 (EASY) Indications Relating to Deposited Microorganism(s) or Other Biological Material (PCT Rule 13bis)	
0-1-1	Prepared using	PCT-EASY Version 2.9 (updated 01.03.2001)
0-2	International Application No.	
0-3	Applicant's or agent's file reference	000505 BT

1	The indications made below relate to the deposited microorganism(s) or other biological material referred to in the description on:	
1-1	page	18
1-2	line	8-13
1-3	Identification of Deposit	
1-3-1	Name of depositary institution	DSMZ-Deutsche Sammlung von Mikroorganismen und Zellkulturen GmbH
1-3-2	Address of depositary institution	Mascheroder Weg 1b, D-38124 Braunschweig, Germany
1-3-3	Date of deposit	16 January 2001 (16.01.2001)
1-3-4	Accession Number	DSMZ 13994
1-4	Additional Indications	NONE
1-5	Designated States for Which Indications are Made	all designated States
1-6	Separate Furnishing of Indications These indications will be submitted to the International Bureau later	NONE

## FOR RECEIVING OFFICE USE ONLY

0-4	This form was received with the international application: (yes or no)	yes
0-4-1	Authorized officer	 B. GATINET (0)70/3402181

## FOR INTERNATIONAL BUREAU USE ONLY

0-5	This form was received by the International Bureau on:	
0-5-1	Authorized officer	

BUDAPEST TREATY ON THE INTERNATIONAL  
RECOGNITION OF THE DEPOSIT OF MICROORGANISMS  
FOR THE PURPOSES OF PATENT PROCEDURE

INTERNATIONAL FORM

Degussa-Hüls AG  
Kantstr. 2  
33790 Halle/Künsebeck

RECEIPT IN THE CASE OF AN ORIGINAL DEPOSIT  
issued pursuant to Rule 7.1 by the  
INTERNATIONAL DEPOSITARY AUTHORITY  
identified at the bottom of this page

<b>I. IDENTIFICATION OF THE MICROORGANISM</b>	
Identification reference given by the DEPOSITOR:  DM1547	Accession number given by the INTERNATIONAL DEPOSITARY AUTHORITY:  DSM 13994
<b>II. SCIENTIFIC DESCRIPTION AND/OR PROPOSED TAXONOMIC DESIGNATION</b>	
The microorganism identified under I. above was accompanied by:  <div style="margin-left: 40px;"> <input type="checkbox"/> a scientific description  <input checked="" type="checkbox"/> a proposed taxonomic designation         </div> (Mark with a cross where applicable).	
<b>III. RECEIPT AND ACCEPTANCE</b>	
This International Depositary Authority accepts the microorganism identified under I. above, which was received by it on 2001-01-16 (Date of the original deposit) <sup>1</sup> .	
<b>IV. RECEIPT OF REQUEST FOR CONVERSION</b>	
The microorganism identified under I above was received by this International Depositary Authority on (date of original deposit) and a request to convert the original deposit to a deposit under the Budapest Treaty was received by it on (date of receipt of request for conversion).	
<b>V. INTERNATIONAL DEPOSITARY AUTHORITY</b>	
Name: DSMZ-DEUTSCHE SAMMLUNG VON MIKROORGANISMEN UND ZELLKULTUREN GmbH  Address: Mascheroder Weg 1b D-38124 Braunschweig	Signature(s) of person(s) having the power to represent the International Depositary Authority or of authorized official(s):  <div style="text-align: center; font-family: cursive; font-size: 1.2em;">U. Wicks</div> Date: 2001-01-18

<sup>1</sup> Where Rule 6.4 (d) applies, such date is the date on which the status of international depositary authority was acquired.


BUDAPEST TREATY ON THE INTERNATIONAL  
RECOGNITION OF THE DEPOSIT OF MICROORGANISMS  
FOR THE PURPOSES OF PATENT PROCEDURE

INTERNATIONAL FORM

Degussa-Hüls AG  
Kantstr. 2  
33790 Halle/Künsebeck

VIABILITY STATEMENT

issued pursuant to Rule 10.2 by the  
INTERNATIONAL DEPOSITARY AUTHORITY  
identified at the bottom of this page

<b>I. DEPOSITOR</b>  Name: Degussa-Hüls AG Kantstr. 2 Address: 33790 Halle/Künsebeck	<b>II. IDENTIFICATION OF THE MICROORGANISM</b>  Accession number given by the INTERNATIONAL DEPOSITARY AUTHORITY: DSM 13994  Date of the deposit or the transfer <sup>1</sup> : 2001-01-16
<b>III. VIABILITY STATEMENT</b>  The viability of the microorganism identified under II above was tested on 2001-01-16 <sup>2</sup> . On that date, the said microorganism was  <input checked="" type="checkbox"/> viable <input type="checkbox"/> no longer viable	
<b>IV. CONDITIONS UNDER WHICH THE VIABILITY TEST HAS BEEN PERFORMED<sup>4</sup></b>  <div style="height: 40px; border: 1px solid black;"></div>	
<b>V. INTERNATIONAL DEPOSITARY AUTHORITY</b>	
Name: DSMZ-DEUTSCHE SAMMLUNG VON MIKROORGANISMEN UND ZELLKULTUREN GmbH  Address: Mascheroder Weg 1b D-38124 Braunschweig	Signature(s) of person(s) having the power to represent the International Depositary Authority or of authorized official(s):    Date: 2001-01-18

- <sup>1</sup> Indicate the date of original deposit or, where a new deposit or a transfer has been made, the most recent relevant date (date of the new deposit or date of the transfer).
- <sup>2</sup> In the cases referred to in Rule 10.2(a) (ii) and (iii), refer to the most recent viability test.
- <sup>3</sup> Mark with a cross the applicable box.
- <sup>4</sup> Fill in if the information has been requested and if the results of the test were negative.

## INTERNATIONAL SEARCH REPORT

Int Application No  
PCT/EP 01/10211

A. CLASSIFICATION OF SUBJECT MATTER  
 IPC 7 C12N15/54 C12N15/74 C12N1/21 C12N9/12 C12P13/04  
 C12P13/08 C12Q1/68

According to International Patent Classification (IPC) or to both national classification and IPC

## B. FIELDS SEARCHED

Minimum documentation searched (classification system followed by classification symbols)  
 IPC 7 C12N C07K

Documentation searched other than minimum documentation to the extent that such documents are included in the fields searched

Electronic data base consulted during the International search (name of data base and, where practical, search terms used)

EPO-Internal, WPI Data, PAJ, BIOSIS, CHEM ABS Data

## C. DOCUMENTS CONSIDERED TO BE RELEVANT

Category *	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
X	DATABASE SWALL 'Online! EBI; 1 October 1996 (1996-10-01) "M. leprae PknB" XP002185694 Acc. No. P54744	1-7,20, 21
P,X	EP 1 108 790 A (KYOWA HAKKO KOGYO KK) 20 June 2001 (2001-06-20) SEQ ID NO:3546 table 1	1-24
A	EP 1 029 919 A (DEGUSSA ;KERNFORSCHUNGSANLAGE JUELICH (DE)) 23 August 2000 (2000-08-23) abstract	8-19

☒ Further documents are listed in the continuation of box C.

☒ Patent family members are listed in annex.

## \* Special categories of cited documents:

- \*A\* document defining the general state of the art which is not considered to be of particular relevance
- \*E\* earlier document but published on or after the International filing date
- \*L\* document which may throw doubts on priority claim(s) or which is cited to establish the publication date of another citation or other special reason (as specified)
- \*O\* document referring to an oral disclosure, use, exhibition or other means
- \*P\* document published prior to the international filing date but later than the priority date claimed

\*T\* later document published after the international filing date or priority date and not in conflict with the application but cited to understand the principle or theory underlying the invention

\*X\* document of particular relevance; the claimed invention cannot be considered novel or cannot be considered to involve an inventive step when the document is taken alone

\*Y\* document of particular relevance; the claimed invention cannot be considered to involve an inventive step when the document is combined with one or more other such documents, such combination being obvious to a person skilled in the art.

\*G\* document member of the same patent family

Date of the actual completion of the international search

14 December 2001

Date of mailing of the international search report

11/01/2002

Name and mailing address of the ISA

European Patent Office, P.B. 5618 Patentlaan 2  
 NL - 2280 HV Rijswijk  
 Tel (+31-70) 340-2040, Tx. 31 651 epo nl,  
 Fax (+31-70) 340-3016

Authorized officer

Mata Vicente, T.

# INTERNATIONAL SEARCH REPORT

Im Application No  
PCT/EP 01/10211

C.(Continuation) DOCUMENTS CONSIDERED TO BE RELEVANT		
Category *	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
A	<p>CREMER J ET AL: "CONTROL OF THE LYSINE BIOSYNTHESIS SEQUENCE IN CORYNEBACTERIUM GLUTAMICUM AS ANALYZED BY OVEREXPRESSION OF THE INDIVIDUAL CORRESPONDING GENES" APPLIED AND ENVIRONMENTAL MICROBIOLOGY, WASHINGTON,DC, US, vol. 57, no. 6, 1 June 1991 (1991-06-01), pages 1746-1752, XP000616281 ISSN: 0099-2240</p>	8-19

# INTERNATIONAL SEARCH REPORT

International Application No  
PCT/EP 01/10211

Patent document cited in search report		Publication date	Patent family member(s)	Publication date
EP 1108790	A	20-06-2001	EP 1108790 A2	20-06-2001
EP 1029919	A	23-08-2000	DE 19907347 A1	24-08-2000
			AU 1760000 A	24-08-2000
			BR 0000897 A	02-05-2001
			CN 1266905 A	20-09-2000
			EP 1029919 A2	23-08-2000
			JP 2000236893 A	05-09-2000
			SK 1942000 A3	12-09-2000